F

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GB

1 Attaching the support handle to the wall

1 Slide the painted plate attachment cover (B) over the handle (A). 2 Fit the wall attachment (C) in Ultraform® to the handle with the special stainless steel screw (D) and tighten with a 6 Allen spanner. 3 Use the wall attachment (C) as a template for tracing the drilling points on the wall. 4 Fix it all to the wall with suitable expansion bolts (see the "anchoring systems" section). 5 Cover the wall attachment (C) with the attachment cover (B).

2 Fixing the tilting support to the wall

1 Use the paper fixing template as a reference for tracing drilling points on the wall. 2 Fix it all to the wall with suitable expansion bolts (see the "anchoring systems" section). It is possible to adjust the resistance to movement of the bar by opening or closing the stainless steel Allen screws (A).

3 Use of the attachment cover with hole

If tightening of the screws within the expansion bolts is problematic, an attachment cover with hole has been supplied for inserting a screwdriver. 1 Slide the painted plate attachment cover with hole (A) over the handle. 2 Fit the wall attachment in Ultraform® to the handle with the special stainless steel screw and tighten with a 6 Allen spanner. 3 Cover the wall attachment with the attachment cover (A). 4 Turn the attachment cover (A) to find the holes for fixing the expansion bolts and mark the drill points on the wall. 5 After inserting the expansion bolts in the wall, turn the attachment cover (A) until it comes into contact with the expansion bolt and screw the screw by inserting a screwdriver in the hole on the attachment cover. 6 Close the hole with the supplied plug (B).

4 Anchoring systems

The wall anchoring system of our handles has been designed to be able to withstand loads to 150 kg. However, based on studies effected with Polytechnic University of Milan, special attention needs to be taken in selecting the anchoring system based on the type of wall in order not to compromise the operation of the device. First of all the selected expansion bolts need to be of a good quality and the wall must have a good resistance. We have provided some installation examples below. **Example 1 - Concrete, cellular concrete or brick.** Use 8 x 80 wood screws and 10 mm universal expansion bolts. **Example 2 - Sheetrock with reinforcement.** Use 3 hollow wall expansion bolts. **Example 3 - Sheetrock with reinforcement and rear support.** To obtain greater wall stability with lightweight material, use a weight bearing reinforcement or add a rear facing. If working on walls which are not compact, it is advisable to use at least one chemical expansion bolt for each wall attachment, if in doubt contact a supplier specialised in anchoring systems.

Intended use. Bathrooms intended for use by the disabled (Italian Ministerial Decree 236/89 and Italian Presidential Decree 503/96) or elderly individuals. Our devices have been tested by Polytechnic University of Milan in accordance with the standard UNI 11064 "Handles: safety requirements and testing methods" and are able to withstand a load up to 150 kg if suitably anchored to walls or weight bearing structures.

Declaration of conformity. Our products are manufactured with materials which are non-toxic, fire proof and non-irritating or sensitising upon contact with the skin. They are in compliance with Community Directive no. 93/42 related to "Medical Devices".

Cleaning and hygiene. Use common detergents and sanitising agents. Do not use powders even if diluted, products with a strong acid or strong oxidising base which could damage the aesthetic or functional characteristics of the product.

Guarantee. The conditions for warranties against manufacturing defects in the country where you purchased this product are valid. Please contact your retailer for additional information.

Warning. Poor maintenance may result in an unstable hold, while poor attachment to the wall may cause the device to detach and the user to fall.

In order to prevent water from accumulating in the seat, please ensure it is dried and returned to the upright position after use.

Strictly observe the safety precautions and operating instructions.

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